**INSTRUCTIONS FOR USE**  
PRESSURE CABLES

-  **TEMPERATURE LIMITS**
-  **KEEP AWAY FROM DIRECT SUNLIGHT**
-  **PROTECT FROM MOISTURE**
-  **LOT CODE**
-  **QUANTITY**
-  **EXPIRY**
-  **DO NOT USE IF PACKAGING IS DAMAGED**
-  **MEDICAL DEVICE**
-  **SEPARATE COLLECTION AS WASTE ELECTRICAL AND ELECTRONIC EQUIPMENT**
-  **ROHS COMPLIANT**
-  **NON-STERILE**
-  **MANUFACTURER**
-  **LATEX FREE**  
Not made with natural rubber latex
-  **PHT FREE**
-  **READ INSTRUCTIONS FOR USE**
-  **REFERENCE**
-  **EC ACCORDING TO EU MDR 2017/745**

**PRESSURE CABLES**

**I. IDENTIFICATION - APPLICATIONS**

The Nissha Medical Technologies Pressure Monitoring Cable plays a critical part in the transmission of the “invasive” blood pressure electric signal. It links the single use pressure sensor (transducer) to the appropriate monitor.

**IMPORTANT**

Use of the pressure cable by unqualified personnel is prohibited.

Also refer to the directions for use of the monitor to which the Nissha Medical Technologies Pressure Monitoring Cable is connected.

**II. MECHANICAL INTEGRITY**

The insulation on a Nissha Medical Technologies Pressure Monitoring Cable and on its conductors is made of polyurethane, giving it excellent tensile and bending strength; most of our connectors are molded and have flexible sleeves to minimize the risk of breakage of the conductors.

**III. SAFETY**

The Nissha Medical Technologies Pressure Monitoring Cable connection system does not include any accessible metallic parts.

The contact pins are gold-plated for better contact and are recessed to protect patients from any contact between these pins and live external metallic parts.

**IV. CLEANING / DISINFECTION / STERILIZATION**

**CLEANING**

As often as necessary, clean the insulating parts of the Nissha Medical Technologies Pressure Monitoring Cable with warm soapy water, then rinse thoroughly.

**DISINFECTION**

Intensive use of the Nissha Medical Technologies Pressure Monitoring Cable requires disinfecting at least once a week in accordance with the following recommendations, taken from the study of the ANIOS laboratory (reference no. 64 16.94/0387).

1. Fill a soaking tank with an 0.5% solution of HEXANIOS G+R solution.
2. Partially immerse the cable (protect the ends to avoid any electrical trouble with the connectors).
3. Allow a contact time of 15 minutes.
4. Rinse the cable.
5. Dry the cable with absorbent tissue.
6. During the immersion, rub the ends of the cable with a LINGET'ANIOS.

**N.B.:** Change the solution every 48 hours.

**STERILIZATION**

A sterilization procedure may have to be performed if the Nissha Medical Technologies Pressure Monitoring Cables are used in a sterile room. It is recommended that the cables be sterilized by the following process after they have been disinfected as described above.

- sterilization time: 22 hours, including 20 hours' exposure to the gas
- sterilizing agent: carbon dioxide + ethylene oxide (80/20%)
- initial vacuum: -70 kPa
- relative humidity: > 60%
- temperature: 50°C
- concentration of ethylene oxide: 530 g/m<sup>3</sup>
- pressure: 120 kPa
- final vacuum: -70 kPa
- rinsing

Investigation protocol 003: sterilization test/pressure cables - AQ/VD/IP/SC/PI - Revision 0 of 3/06/98  
Never sterilize the pressure cables with steam or boiling water.

**V. INSTALLATION INSTRUCTIONS**

Check the reference of the Nissha Medical Technologies Pressure Monitoring Cable, marked on the packaging or on the cable itself, to make sure that it is compatible with the monitor. To disconnect the cable from the monitor, pull on the connector, not the cable itself.

For the connections of the cable to the monitor and to the pressure sensor, refer to their respective directions for use. When the Nissha Medical Technologies Pressure Monitoring Cable is not in use, it is important to make sure no liquid gets into the connectors. Moisture in the connector can cause the cable to malfunction or alter the pressure reading.

**VI. WARRANTY**

Unused cables still in their original packaging are guaranteed by Nissha Medical Technologies for one year counting from the date of manufacture, unless damaged by mishandling.

Nissha Medical Technologies cannot guarantee the cable against damage during use.

Nissha Medical Technologies guarantees the compliance of the cable with the specifications of the safety and performance standards applicable to it and currently in force.

**VII. USE - MAINTENANCE**

A reduced operating test is necessary before each use.

**CAUTION:** to avoid potential errors due to a faulty cable, Nissha Medical Technologies recommends the function of the cable be verified prior to actual use. This can be accomplished by using the BD XCALIBERTM as a known pressure source, or equivalent methods.

**CAUTION:** checking the vent tubes is necessary before each use. Pending use, pressure cables must be kept in their original packaging and stored under the following conditions:

- Ambient temperature between +10°C and +40°C

Depending on the number of disinfections/sterilizations:

- Visually check the integrity of the cable's insulating sheath.
- Check the electrical continuity of the conductors of the cable.
- Check that the marking is legible.
- Check the mechanical integrity of the connectors.

If the appearance of one or more parts of the cable is doubtful, the cable must be replaced. Nissha Medical Technologies cannot be held liable for any failure to observe these rules.

Nissha Medical Technologies is not liable for trouble resulting from failure to follow the rules for installing and using the Nissha Medical Technologies Pressure Monitoring Cable described in these directions for use.

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 <b>TEMPERATURE LIMITS</b>	 <b>ROHS COMPLIANT</b>
 <b>KEEP AWAY FROM DIRECT SUNLIGHT</b>	 <b>NON-STERILE</b>
 <b>PROTECT FROM MOISTURE</b>	 <b>MANUFACTURER</b>
 <b>LOT CODE</b>	 <b>LATEX FREE</b>
 <b>QUANTITY</b>	 <b>PHT FREE</b>
 <b>EXPIRY</b>	 <b>READ INSTRUCTIONS FOR USE</b>
 <b>DO NOT USE IF PACKAGING IS DAMAGED</b>	 <b>REFERENCE</b>
 <b>MEDICAL DEVICE</b>	 <b>EC ACCORDING TO EU MDR 2017/745</b>
 <b>SEPARATE COLLECTION AS WASTE ELECTRICAL AND ELECTRONIC EQUIPMENT</b>	

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