



NeuroPlus™

Sleep Study Electrodes

FEATURES: A10043



SUPERIOR ADHESION

Small, aggressive adhesive electrode designed for easy placement and superior adhesion.



SOLID GEL

Specially formulated solid gel ensures full sensor contact and excellent readings.

A10040 | A10041 | A10042:



GENTLE ADHESIVE

Designed with a gentle adhesive gel for use in sensitive areas.



REPOSITIONABLE

Specially formulated adhesive gel is repositionable to ensure proper placement.

For use in sleep studies, EMG, NCS and EP

PACKAGING

RESEALABLE BULK POUCH

Ziplock pouch contains 30 or more electrodes

CASE

FDA / GS1 Compliant Labels



HYGIENIC, DISPOSABLE SINGLE-USE PACKS

Reduces risk of cross contamination of unused electrodes

INNER BOX

Single use packs are packaged in inner boxes within the case

1.800.669.6905 | HS.NisshaMedical.com

Patient Monitoring & Medical Surgical Consumables

Latex-free Solid Gel






NeuroPlus™




Repositionable, Disposable Electrodes

NeuroPlus™ electrodes are smaller in size, allowing for easy placement on the face or other difficult locations. Available in both aggressive and gentle adhesives, NeuroPlus™ features a cloth base or foam option for added comfort to the patient.




A10040 - 0.875" / 22.225mm - SQUARE CLOTH ELECTRODE - REPOSITIONABLE

	Product # A10040-60	Per Strip 3	Per Pouch 60	Per Box n/a	Per Case 600	Solid Gel 	Duration of Use  48 Hours
---	-------------------------------	-----------------------	------------------------	-----------------------	------------------------	---	---




A10041 - 0.875" x 1.3125" / 22.225mm x 33.34mm - RECT. CLOTH ELECTRODE - REPOSITIONABLE

	Product # A10041-60	Per Strip 3	Per Pouch 60	Per Box n/a	Per Case 600	Solid Gel 	Duration of Use  48 Hours
--	-------------------------------	-----------------------	------------------------	-----------------------	------------------------	---	---




A10042 - 1.375" / 34.925mm - SQUARE CLOTH ELECTRODE - REPOSITIONABLE

	Product # A10042-60	Per Strip 3	Per Pouch 60	Per Box n/a	Per Case 600	Solid Gel 	Duration of Use  48 Hours
---	-------------------------------	-----------------------	------------------------	-----------------------	------------------------	---	---

A10043 - 1" / 25.4mm - SQUARE CLOTH ELECTRODE

	Product # A10043	Per Strip 3	Per Pouch 30	Per Box n/a	Per Case 300	Solid Gel 	Duration of Use  48 Hours
---	----------------------------	-----------------------	------------------------	-----------------------	------------------------	---	---

A7-5 - 0.875" / 22.225mm - TEARDROP FOAM ELECTRODE

	Product # A7-5	Per Strip 5	Per Pouch 5	Per Box 75	Per Case 750	Solid Gel 	Duration of Use  72 Hours
---	--------------------------	-----------------------	-----------------------	----------------------	------------------------	---	---