



# Performance Plus™ Electrodes

Electrodes for Stress and Holter Testing

## FEATURES



### HELPS REDUCE ARTIFACTS

Aggressive adhesive keeps electrodes in place on patient - reducing motion artifact and repeat testing.



### HIGH CHLORIDE WET GEL

High chloride wet gel hydrates the skin and produces fast and clear traces.



### Ag/AgCL SENSOR

Silver/Silver Chloride sensor provides accurate trace results.



### HYPOALLERGENIC

Hypoallergenic adhesive, not made with natural rubber latex. Free from common allergens.

In critical diagnostic situations, choose Vermed PerformancePlus™ Electrodes.

## PACKAGING

### RESEALABLE BULK POUCH

Ziplock pouch contains 30 or more electrodes

### CASE

FDA / GS1 Compliant Labels



### HYGIENIC, DISPOSABLE SINGLE-USE PACKS

Reduces risk of cross contamination of unused electrodes

### INNER BOX

Single use packs are packaged in inner boxes within the case

1.800.669.6905 | [HS.NisshaMedical.com](http://HS.NisshaMedical.com)

Patient Monitoring & Medical Surgical Consumables

Latex-free




Wet Gel




# PerformancePlus™ Electrodes

Electrodes for Stress and Holter Testing


## A10005 - 1.5"/38.1mm - TEARDROP FOAM ELECTRODE

	Product #	Per Strip	Per Pouch	Per Box	Per Case	Wet Gel	Duration of Use
	A10005	5	5	100	600		24 Hours
A10005NB	5	5	n/a	600			
A10005-1-60	1	60	n/a	600			
A10005-4NB	4	4	n/a	600			
A10005-7NB	7	7	n/a	700			
A10005-60	5	60	n/a	600			


## A10006-5 - 1.875"/47.625mm - TEARDROP FOAM ELECTRODE

	Product #	Per Strip	Per Pouch	Per Box	Per Case	Wet Gel	Duration of Use
	A10006-5	5	5	100	600		24 Hours
A10006-1-60S	1	60	n/a	600			


## A10006-5T - 2.25"/57.15mm - ROUND TAPE ELECTRODE

	Product #	Per Strip	Per Pouch	Per Box	Per Case	Wet Gel	Duration of Use
	A10006-5T	5	5	100	600		24 Hours
A1006-1-60TS	1	60	n/a	600			


## A10007 - 2"/50.8mm - ROUND CLOTH ELECTRODE

	Product #	Per Strip	Per Pouch	Per Box	Per Case	Wet Gel	Duration of Use
	A10007-5	5	5	100	600		24 Hours

## A10012 - 1.875" x 1.25" / 47.625mm x 31.75mm - RECT. FOAM ELECTRODE

	Product #	Per Strip	Per Pouch	Per Box	Per Case	Wet Gel	Duration of Use
	A10012-5S	5	5	100	600		24 Hours
A10012-4-60S	4	60	n/a	600			
A10012-5S	5	5	100	600			
A10012-6S	6	6	n/a	600			
A10012-60S	5	60	n/a	600			

## A10101-3 - 1.875" x 1.25" / 47.625mm x 31.75mm - RECT. FOAM ELECTRODE

	Product #	Per Strip	Per Pouch	Per Box	Per Case	Wet Gel	Duration of Use
	A10101-10	5	10	100	1000		24 Hours